



Biomedical Systems

Microscope Slides




**SLIDES OPTIMIZED FOR
AEROSPRAY® SLIDE STAINERS,
CYTOCENTRIFUGES, AND
CELLAVISION® DIGITAL
CELL ANALYZERS**


MICROSCOPE BLANK SLIDES

ELITechGroup microscope slides are optimized for Aerospray® Slide Stainers, Cytocentrifuges, and CellaVision® Digital Cell Analyzers.

- Made of white soda lime glass enabling high quality smearing
- Single use



Labeling surface area compatible with pencils, StatPens, and barcodes.



Clipped corner microscope slides (45°) compatible with CellaVision analyzers.

SPECIFICATIONS

| | |
|------------|--|
| Dimensions | 75.0 mm x 25.0 mm (+/- 0.2 mm) |
| Thickness | 1.0 mm thick with a tolerance of +/- 0.10 mm |
| Material | White Soda Lime Glass |
| Quantity | 72 Slides |

ORDERING INFORMATION

| | |
|--------|---|
| SS-280 | Clipped Corner Glass Microscope Slides, Frosted |
| SS-287 | Glass Microscope Slides, Frosted |



WORLDWIDE OFFICES

Please contact your sales representative for terms, conditions and product availability in your country.

Headquarters

Australia
Belgium & Luxembourg
Brazil
France
Italy
New Zealand

+33 1 41 45 07 10

T: +61 1800 815 098
T: +32 9 282 05 31
T: +55 27 3025 1415
T: +33 4 83 36 10 82
T: +39 02 48 40 35 42
T: +64 800 555 611

Germany
Canada
Serbia
The Netherlands
UK
United States

T: +49 30 7809 588-0
T: +1 888 662 4925
T: +381 11 2467119
T: +31 313 430 500
T: +44 1442 869320
T: +1 435 752 6011



AerosprayStaining.com
aerospray@elitechgroup.com