



SEE WHAT YOU NEED TO SEE

## Aerospray<sup>®</sup> Comparison Chart

---

The only slide stainer with built-in cyto centrifugation and programmable staining capabilities.

Innovation with Integrity

# WHICH AEROSPRAY® SLIDE STAINER IS RIGHT FOR YOUR LAB?



Aerospray®  
Cytology



Aerospray®  
Hematology Stat



Aerospray®  
Hematology Pro



Aerospray®  
Gram



Aerospray®  
TB

	MODEL 7522	MODEL 7122	MODEL 7152	MODEL 7322	MODEL 7722
<b>DISCIPLINE</b>	Cytology	Hematology		Microbiology	
<b>STAINING METHOD</b>	Papanicolaou	Romanowsky	Wright-Giemsa, May Grunwald, Rapid Romanowsky	Gram differentiation	Fluorescence or Carbol Fuchsin
<b>NOZZLE CONFIGURATION</b>	2 Nozzles	4 Nozzles	2 Nozzles	6 Nozzles	
<b>FIXATIVE</b>	Ethanol/Carbowax™	Methanol or Aerofix™ fixative		Methanol	Mycohold™
<b>CAROUSEL CAPACITY</b>	1-12 or 1-30 slides	1-12 slides	1-12 or 1-30 slides	1-12 or 1-30 slides	
<b>THROUGHPUT</b>	~200 slides/hour	~350 slides/hour	~250 slides/hour	~300 slides/hour	~300 slides/hour
<b>RUN TIME</b>	8-15 minutes	2-4 minutes	7-12 minutes	6-8 minutes	6-8 minutes
<b>STAINS</b>	EA-50, Orange G, Bluing Agent, Hematoxylin I and II, and Acid Alcohol	Aqueous-based stains: Thiazin, Eosin, and Rinse; Basofix stains available for improved visualization of Basophil granules		Crystal Violet, Iodine, and Decolorizer with either Safranin or Fuchsin	Fluorescence and modified Ziehl-Neelsen primary, decolorizer, and counterstains

## Aerospray Benefits:

- Eliminate risk of cross-contamination
- Require minimal maintenance
- Programmable to suit your staining preferences
- Reduce labor and stain cost per slide
- Facilitate compliance and lab standardization
- Dual-purpose slide stainer and cytocentrifuge

[www.Aerospraystaining.com](http://www.Aerospraystaining.com)

